

Seeding Cartridge Holders for Model 44D Well Chamber

44D-4 Mick Magazine

The 44D-4 source holder for the Mick Magazine will position iodine seeds securely at the sweet spot of the CNMC Model 44D well chamber allowing batch assay. Cartridge placement is quick to minimize exposure to hands. Simply place the loaded cartridge through the hole provided with the end holding the seeds toward the bottom and insert the source holder into the well chamber. The 44D-4 source holder lid includes a center hole with a nylon set screw and is designed to also accept the 44D-1S acrylic tube for seed sources or the 44D-1L acrylic tube for 137 Cs tube sources. This additional feature provides versatility and economy.

44D-6 Índigo EXPRESS

The 44D-6 source holder for the Índigo EXPRESS Seeding Cartridge positions palladium seeds securely at the sweet spot of the CNMC Model 44D well chamber. Cartridge placement is quick to minimize exposure to hands. Simply place the loaded cartridge through the slot provided with the end holding the seeds toward the bottom and insert the source holder into the well chamber. The 44D-6 source holder lid includes a center hole with a nylon set screw and is designed to also accept the 44D-1S acrylic tube for seed sources or the 44D-1L acrylic tube for 137 Cs tube sources. This additional feature provides versatility and economy.

MODEL 44D-4



MODEL 44D-6

